

Goat anti-MYO Polyclonal Antibody

Basic information:

Catalog No.:	UDA9017-9G	Source:	Goat
Clonality:	Polyclonal	Isotype:	IgG
Purity:	>95% as determined by SDS-PAGE.		

Useful Information:

Applications:	LETIA
Specificity:	Recognizes native Human MYO.
Format:	Liquid
Purification:	Affinity purified.
Immunogen:	native Human MYO.
Presentation:	PBS,pH 7.4
Description:	Myoglobin is a small protein containing heme, which mainly exists in striated muscle. It is responsible for oxygen storage and distribution in cardiac muscle cells. MYO, cTnI and CK-MB are three rapid markers of myocardial infarction
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	Test systems and platform change may cause different results, please contact us if you have any questions. This product is sold for research use only. Standard Laboratory Practices should be followed when handling this material.