Goat anti-ApoA1 Immune Serum

Racin	into	mation:
Dasic		mation.

Dasic Information.					
Catalog No.:	UDA9039-9G	Host:	Goat		
Clonality:	Polyclonal				
Useful Information:					
Applications:	TIA				
Specificity:	Recognizes Human ApoA1.				
Format:	Liquid				
Immunogen:	Native Human ApoA1				
Presentation:	Antisera against whole serums are obtained by immunizing host animals with whole serum. This product is reddish transparent liquid, containing 0.03% Proclin300.				
Description:	Apoa-I is the most important structural component of high density Lipoprotein (HDL), accounting for about 70% of its protein content.It plays a key role in reverse cholesterol transport (RCT), a process that transports cholesterol from peripheral tissues to liver metabolism. Apoa-I, a member of the apolipoprotein superfamily, is a cofactor of cholesterol lipoyltransferase (LCAT)				
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.				
Note:	Test systems and platform cha us if youhave any questions. This product is sold for resear be followed when handling th	ch use only. Sta			