

PCT Recombination Protein

Basic information:

Catalog No.:	UDP9025-1	Source:	Human
Expression Host:	<i>Escherichia coli</i>	Type:	Recombinant protein
Purity:	>90 % as determined by SDS-PAGE		

Useful Information:

Construction:

A DNA sequence encoding the PCT Protein was expressed with His-tag in C terminus.

Molecular Mass:

PCT Recombination Protein consists of 116 amino acids.

Formulation:

Liquidin PBS, PH7.4, containing 20% glycerol

Biological Activity

N/A

Endotoxin:

N/A

Storage:

Recombinant proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

Reconstitution:

According to the application.

Description:

Procalcitonin (PCT) is a precursor peptide of the hormone calcitonin (CT). The content of PCT in the serum of normal people is very low. After being infected by microorganisms, calC-I gene is induced to express and promote the secretion of PCT, reaching the peak at 8 hours. PCT has been recognized as a marker of systemic inflammatory response syndrome to guide treatment, reduce antibiotic use, diagnose sepsis, and improve long-term outcomes.

Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes. Not for human use Made in China.