

NGAL Recombination Protein

Basic information:

Catalog No.:	UDP9029-1	Source:	Human
Expression Host:	<i>Escherichia coli</i>	Type:	Recombinant protein
Purity:	> 90 % as determined by SDS-PAGE		

Useful Information:

Construction:

A DNA sequence encoding the NGAL Protein was expressed with His-tag in N and Cterminus.

Molecular Mass:

The protein of NGAL consists of 178 amino acids.

Formulation:

Liquid in PBS, pH7.4

Biological Activity

N/A

Endotoxin:

N/A

Storage:

Recombinant proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

Reconstitution:

According to the application.

Description:

Neutrophil gelatinase-associated Lipocalin (NGAL) is a novel and highly effective clinical marker for early AKI renal injury.

Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes. Not for human use. Made in China.