

## AMH Recombination Protein

### Basic information:

<b>Catalog No.:</b>	UDP9034-2	<b>Source:</b>	Human
<b>Expression Host:</b>	HEK293	<b>Type:</b>	Recombinant protein
<b>Purity:</b>	> 90 % as determined by SDS-PAGE		

### Useful Information:

#### Construction:

A DNA sequence encoding the AMH Protein was expressed with His-tag in C terminus.

#### Molecular Mass:

The protein of AMH consists of 535 amino acids.

#### Formulation:

Liquidin PBS, pH7.4

#### Biological Activity

N/A

#### Endotoxin:

N/A

#### Storage:

Recombinant proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

#### Reconstitution:

According to the application.

#### Description:

AMH provides a reliable and rapid assessment of ovarian reserve function, the number of eggs likely to form in a woman's ovaries, and thus her fertility and the onset of menopause.

#### Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes. Not for human use. Made in China.