# **ApoB Recombination Protein**

**Basic information:** 

Catalog No.: UDP9097-2 Source: Human

**Expression Host:** HEK293 **Type:** Recombinant protein

**Purity:** >90 % as determined by SDS-PAGE

# **Useful Information:**

#### **Construction:**

A DNA sequence encoding the ApoB Protein was expressed with His-tag in C terminus.

#### **Molecular Mass:**

ApoB Recombination Protein consists of 298 amino acids.

#### Formulation:

Liquidin PBS,PH7.4

# **Biological Activity**

N/A

#### **Endotoxin:**

N/A

# Storage:

Recombinant proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

# **Reconstitution:**

According to the application.

# **Description:**

ApoB protein is the most important apolipoprotein on chylomicron and LOW density lipoprotein, which can form low density lipoprotein cholesterol and remove the residue of lipid metabolism through endocytosis mediated by low density lipoprotein receptor (LDLR) on the surface of liver cells.

### Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes. Notfor human use Made in China.