

Native Human Urine mALB Protein

Basic information:

Catalog No.:	UDP9107-9	Species:	Human
Expression Host:	<i>Native</i>	Type:	Full length protein

Useful Information:

Purity:

> 90 % as determined by SDS-PAGE

Formulation:

Liquidin PBS pH7.4, containing 10% glycerin and 0.05% Proclin 300.

Biological Activity

N/A

Endotoxin:

N/A

Storage:

Proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

Reconstitution:

According to the application.

Description:

The determination of microalbumin in urine reflects early nephropathy and renal injury.

Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes.
Not for human use. Made in China.