Goat anti-Human IgG Immune Serum

Basic information:

Catalog No.: UB0009-9G Host: Goat

Clonality: Polyclonal

Useful Information:

Applications: TIA

Specificity: Recognizes Human IgG.

Format: Liquid

Immunogen: Native Human IgG

Presentation: Antisera against whole serums are obtained by immunizing host animals with

whole serum. This product is reddish transparent liquid, containing 0.03%

Proclin300.

Description: IgG is a monomeric immunoglobulin, built of two heavy chains gamma and two

light chains. Each molecule has two antigen binding sites. This is the most abundant immunoglobulin and is approximately equally distributed in blood and in tissue liquids, constituting 75% of serum immunoglobulins in humans.

There are 4 subclasses: IgG1 (66%), IgG2 (23%), IgG3 (7%) and IgG4 (4%).

Storage: Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.

Note: Test systems and platform change may cause different results, please contact

us if youhave any questions.

This product is sold for research use only. Standard Laboratory Practices should

be followed when handling this material.