Goat anti-Human IgM Immune Serum

Basic information:

Catalog No.: UB0010-9G Host: Goat

Clonality: Polyclonal

Useful Information:

Applications: TIA

Specificity: Recognizes Human IgM.

Format: Liquid

Immunogen: Native Human IgM

Presentation: Antisera against whole serums are obtained by immunizing host animals with

whole serum. This product is reddish transparent liquid, containing 0.03%

Proclin300.

Description: IgM normally constitutes about 10% of serum immunoglobulins. IgM antibody

is prominent in early immune responses to most antigens and predominates in

certain antibody responses such as 'natural' blood group antibodies.

Storage: Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.

Note: Test systems and platform change may cause different results, please contact

us if youhave any questions.

This product is sold for research use only. Standard Laboratory Practices should

be followed when handling this material.