

Native Human IgM protein

Basic information:

Catalog No.:	UB0010	Species:	Human
Expression Host:	<i>Native</i>	Type:	Full length protein

Useful Information:

Purity:

>80% as determined by SDS-PAGE

Formulation:

Liquidin PBS pH7.4, containing 10% glycerin and 0.05% Proclin 300.

Biological Activity

N/A

Endotoxin:

N/A

Storage:

Proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

Reconstitution:

According to the application.

Description:

IgM normally constitutes about 10% of serum immunoglobulins. IgM antibody is prominent in early immune responses to most antigens and predominates in certain antibody responses such as 'natural' blood group antibodies.

Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes.

Not for human use.

Made in China.