## **Goat anti-ApoB Immune Serum**

**Basic information:** 

Catalog No.: UDA9040-9G Host: Goat

**Clonality:** Polyclonal

**Useful Information:** 

Applications: TIA

**Specificity:** Recognizes Human ApoB.

Format: Liquid

Immunogen: Native Human ApoB

Presentation: Antisera against whole serums are obtained by immunizing host animals with

whole serum. This product is reddish transparent liquid, containing 0.03%

Proclin300.

**Description:** ApoB protein is the most important apolipoprotein on chylomicron and LOW

density lipoprotein, which can form low density lipoprotein cholesterol and remove the residue of lipid metabolism through endocytosis mediated by low

density lipoprotein receptor (LDLR) on the surface of liver cells.

**Storage:** Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.

**Note:** Test systems and platform change may cause different results, please contact

us if youhave any questions.

This product is sold for research use only. Standard Laboratory Practices should

be followed when handling this material.