Goat anti-Human mALB Immune Serum

Basic information:

Catalog No.: UDA9107-9G Host: Goat

Clonality: Polyclonal

Useful Information:

Applications: TIA

Specificity: Recognizes Human mALB.

Format: Liquid

Immunogen: Native Human mALB

Presentation: Antisera against whole serums are obtained by immunizing host animals with

whole serum. This product is reddish transparent liquid, containing 0.03%

Proclin300.

Description: The determination of microalbumin in urine reflects early nephropathy and

renal injury.

Storage: Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.

Note: Test systems and platform change may cause different results, please contact

us if youhave any questions.

This product is sold for research use only. Standard Laboratory Practices should

be followed when handling this material.