

Goat anti-Human FN Immune Serum

Basic information:

Catalog No.:	UDA9109-9G	Host:	Goat
Clonality:	Polyclonal		

Useful Information:

Applications:	TIA
Specificity:	Recognizes Human FN.
Format:	Liquid
Immunogen:	Native Human FN
Presentation:	Antisera against whole serums are obtained by immunizing host animals with whole serum. This product is reddish transparent liquid, containing 0.03% Proclin300.
Description:	Fibronectin(FN) is a high molecular weight non-collagenous glycoprotein with 4.5-9.5% sugar content. It is mainly produced by liver and vascular endothelial cells and widely exists in plasma, various cell surfaces and cell matrix. FN plays a key role in regulating cell adhesion, migration, growth, proliferation and differentiation, and is involved in wound healing, embryo development, immune regulation and other processes, with a variety of biological activities and biological functions of cell regeneration and repair.
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.
Note:	Test systems and platform change may cause different results, please contact us if you have any questions. This product is sold for research use only. Standard Laboratory Practices should be followed when handling this material.