## Goat anti-Human PA Immune Serum

**Basic information:** 

Catalog No.: UDA9114-9G Host: Goat

**Clonality:** Polyclonal

**Useful Information:** 

Applications: TIA

**Specificity:** Recognizes Human Prealbumin protein.

Format: Liquid

Immunogen: Native Human Prealbumin protein.

Presentation: Antisera against whole serums are obtained by immunizing host animals with

whole serum. This product is reddish transparent liquid, containing 0.03%

Proclin300.

**Description:** Prealbumin is called thyroxine protein, with a molecular weight of 54,000.It is

synthesized by liver cells. In electrophoretic separation, it is often shown in front of albumin and has a short half-life of only about 1.9 days. Prealbumin has a higher sensitivity than albumin and transferrin when protein is

malnourished and liver dysfunction occurs.

**Storage:** Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.

**Note:** Test systems and platform change may cause different results, please contact

us if youhave any questions.

This product is sold for research use only. Standard Laboratory Practices should

be followed when handling this material.