

ST2s Recombination Protein

Basic information:

Catalog No.:	UDP9047-2	Source:	Human
Expression Host:	HEK293 Cells	Type:	Recombinant protein
Purity:	> 90 % as determined by SDS-PAGE		

Useful Information:

Construction:

A DNA sequence encoding the ST2 Protein was expressed with His-tag in C terminus.

Molecular Mass:

The protein of ST2 consists of 310 amino acids.

Formulation:

Liquidin PBS, pH7.4.

Biological Activity

N/A

Endotoxin:

N/A

Storage:

Recombinant proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

Reconstitution:

According to the application.

Description:

ST2 is a member of the interleukin-1 receptor/Toll-like receptor superfamily. In 2013, THE ACCF/AHA guidelines included sST2 as a marker of myocardial fibrosis that can predict admission and death in HF patients.

Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes. Not for human use. Made in China.