

AQP4 Recombination Protein

Basic information:

Catalog No.:	UDP9104-1	Source:	Human
Expression Host:	<i>Escherichia coli</i>	Type:	Recombinant protein
Purity:	> 90 % as determined by SDS-PAGE		

Useful Information:

Construction:

A DNA sequence encoding the human AQP4 protein was expressed with Trx and His-tag in N terminus.

Molecular Mass:

The protein consists of 241 amino acids.

Formulation:

Liquid in PBS, pH7.4.

Biological Activity

N/A

Endotoxin:

N/A

Storage:

Recombinant proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

Reconstitution:

According to the application.

Description:

Aquaporin-4 (AQP4) is a water channel expressed in astrocytes throughout the central nervous system, as well as in epithelial cells in various peripheral organs. AQP4 is involved in brain water balance, neuroexcitation, astrocyte migration, and neuroinflammation and is the target of pathogenic autoantibodies in neuromyelitis optica. Two AQP4 isoforms produced by alternative splicing, M1 and M23 AQP4, form heterotetramers that assemble in cell plasma membranes in supramolecular aggregates called orthogonal arrays of particles (OAPs).

Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes. Not for human use. Made in China.