

## Prealbumin Native protein

### Basic information:

<b>Catalog No.:</b>	UDP9114-9	<b>Species:</b>	Human
<b>Expression Host:</b>	<i>Native</i>	<b>Type:</b>	Full length protein

### Useful Information:

#### Purity:

> 90 % as determined by SDS-PAGE

#### Formulation:

Liquidin PBS pH7.4, containing 5% glycerin and 0.01% Sodium azide.

#### Biological Activity

N/A

#### Endotoxin:

N/A

#### Storage:

Proteins are provided as frozen liquid which are shipped with dry ice. Bulk packages can be provided as lyophilized powder which shipped with blue ice.

#### Reconstitution:

According to the application.

#### Description:

Prealbumin is called thyroxine protein, with a molecular weight of 54,000. It is synthesized by liver cells. In electrophoretic separation, it is often shown in front of albumin and has a short half-life of only about 1.9 days. Prealbumin has a higher sensitivity than albumin and transferrin when protein is malnourished and liver dysfunction occurs

#### Note:

This material is offered by Gene Universal for research, laboratory or further evaluation purposes.  
Not for human use. Made in China.