Goat anti-Chicken IgY polyclonal antibody

Basic information:

Catalog No.:UB0005-1Source:GoatClonality:PolyclonalIsotype:IgG

Purity: >95% as determined by SDS-PAGE.

Useful Information:

Applications: ELISA,Immunomicroscopy,WB, Dot blot, ICA

Specificity: Recognizes native Chicken IgY.

Format: Liquid

Purification: Affinity purified.

Immunogen: Native Chicken IgY.

Presentation: PBS,pH 7.4.

Description: IgY has atetrameric structure consisting of two identical light chains (Mr = 24

kDa) and two identical heavychains (Mr = 73 kDa). Goat anti-Chicken IgY

polyclonal antibody can recognizeslight chains and heavychains.

Storage: Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.

Note: Test systems and platform change may cause different results, please contact

us if youhave any questions.

This product is sold for research use only. Standard Laboratory Practices should

be followed when handling this material.