## Rabbit anti-Chicken IgY polyclonal antibody

**Basic information:** 

Catalog No.:UB0005-2Source:RabbitClonality:PolyclonalIsotype:IgG

**Purity:** >95% as determined by SDS-PAGE.

**Useful Information:** 

**Applications:** ELISA, Immunomicroscopy, WB, Dot blot, IHC/ICC

**Specificity:** Recognizes native Chicken IgY.

Format: Liquid

Purification: Affinity purified.

Immunogen: Native Chicken IgY.

**Presentation:** PBS PH7.4

**Description:** Secreted immunoglobulins have two separable functions: one is to bind

specifically to pathogenicmolecules, the other is to recruit cells and molecules to destroy the opsonized pathogen. IgY has atetrameric structure consisting of two identical light chains (Mr=24 kDa) and two identical heavychains (Mr=73

kDa).

**Storage:** Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.

**Note:** Test systems and platform change may cause different results, please contact

us if youhave any questions.

This product is sold for research use only. Standard Laboratory Practices should

be followed when handling this material.