Goat anti-Mouse IgG Polyclonal Antibody

Basic information:			
Catalog No.:	UB0012-1	Source:	Goat
Clonality:	Polyclonal	Isotype:	lgG
Purity:	>95% as determined by SDS-PAGE.		
Useful Information:			
Applications:	ELISA, Immunomicroscopy, WB, Dot blot, ICA		
Specificity:	Recognizes native Mouse IgG.		
Format:	Liquid		
Purification:	Affinity purified.		
Immunogen:	Native Mouse IgG.		
Presentation:	PBS,pH 7.4, containing 15% glycerin.		
Description:	ion: Goat anti Mouse IgG antibody recognises the heavy and light chains of murine		
	IgG1, IgG2a, IgG2b and IgG3 as demonstrated by ELISA and Flow Cytometry.		
This antibody may also recognises the light chains			of murine IgA and IgM.
Storage:	Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.		