ST2s (7B10) Monoclonal Antibody

Basic information:

Catalog No.:UDA9047-1Source:MouseClonality:MonoclonalIsotype:IgG2b

Purity: >95% as determined by SDS-PAGE

Useful Information:

Applications: ELISA, CLIA

Usage: Capture or Detection

Specificity: Recognizes ST2s Recombination Protein

Format: Liquid

Purification: Protein A/G affinity purification from mouse ascites.

Immunogen: Recombinant protein of ST2s

Presentation: PBS, pH 7.4

Description: ST2 is a member of the interleukin-1 receptor/Toll-like receptor superfamily.In

2013, THE ACCF/AHA guidelines included sST2 as a marker of myocardial

fibrosis that can predict admission and death in HF patients.

Storage: Store at 4°C short term and -20°C long term. Avoid freeze-thaw cycles.

Note: Test systems and platform change may cause different results, please contact

us if youhave any questions.

This product is sold for research use only. Standard Laboratory Practices should

be followed when handling this material.